



Can I move my PATGuard 2 licence to another PC?

It is possible to
select Trans

1. The Trans

activated.



2. Another computer with PATGuard 2 already installed will now be able to run the program. To return to the original computer, select Transfer... to deactivate the licence. At the original computer select Restore... to re-activate the licence and run



The above

process requires internet connection.

If you have lost your software disk please reinstall using the below link.

www.seaward.com/software/PATGuard_2_1.exe

You will need your serial number to activate the software on the new PC. This will have been on the inside cover of the DVD style software case. If you no longer have this then please follow the link below.

[I've lost my PATGuard Serial Number and/or Unlock Code. What do I do?](#)

If you require more help, please contact us at

<https://www.seaward.com/gb/enquire/>.